




Kyphoplasty System



Concentration Makes Profession

Kyphoplasty Tool Kit - 8G



KT-01-01
Puncture device




KT-02-01
Puncture expander module




KT-03-01
Expander device



KT-04-01
Bone cement filling device



KT-05-01
Bone drill




KT-06-01
Guide wire

Intended use


- This device is mainly used for bone percutaneous intervention and establishing working channel by medical institutions.

Kyphoplasty Tool Kit - 11G




KT-07-01

Puncture expander module



KT-08-01

Bone cement filling device



KT-09-01

Bone drill

Features

- Minimally invasive surgery can reduce tissue injury and blood loss.
- The operation is simple, fast, and the operation time greatly shortened.
- Local anesthesia can be adopted to reduce the risk of surgery and postoperative complications.
- The product can immediately relieve postoperative pain, effectively restore vertebral body height and correct spine kyphosis.

Kyphoplasty Balloon Catheter

Intended use

- This device is mainly used in percutaneous kyphoplasty(pkp) operation to dilate vertebral body and form a cavum which is for injecting bone cement to recover and stabilize vertebral body.

SizeType	Actual Photo
8G	
11G	

Specifications

Model	Distance of Two Radiopaque Markers	Channel ID	Overall Length	Maximum Volume	Constrained Burst Pressure	SizeType
KB0210	10	≥3.65mm	315mm	4cc	≥400 PSI	8G
KB0115	15	≥3.65mm	315mm	4cc	≥400 PSI	8G
KB0120	20	≥3.65mm	315mm	6cc	≥400 PSI	8G
KB0210S1	10	≥3.10mm	280mm	3cc	≥400 PSI	11G
KB0115S1	15	≥3.10mm	280mm	4cc	≥400 PSI	11G
KB0120S1	20	≥3.10mm	280mm	6cc	≥400 PSI	11G

- Customization is also acceptable.


Kyphoplasty tool kit components list - 8G

Serial number	Name	Specification	KT-00-01	KT-00-02	KT-00-03	KT-00-04	KT-00-16	KT-00-17
1	Puncture device	KT-01-01	2	2	1	1	/	/
2	Puncture expander module	KT-02-01	1	/	1	/	2	1
3	Expander device	KT-03-01	2	2	1	1	/	/
4	Bone cement filling device	KT-04-01	6	6	3	3	6	3
5	Bone drill	KT-05-01	1	1	1	1	1	1
6	Guide wire	KT-06-01	2	2	1	1	/	/

Kyphoplasty tool kit components list - 11G

Serial number	Name	Specification	KT-00-05	KT-00-06
1	Puncture expander module	KT-07-01	2	1
2	Bone cement filling device	KT-08-01	6	3
3	Bone drill	KT-09-01	1	1

Balloon Inflator



Intended use

- This device is mainly used for hydraulic pressure on balloons to achieve the function of balloon dilation.

Features

- Optimal head fluid channel design makes it convenient for venting.
- O-type plunger design makes for low frictional resistance and can pressure easily.
- Screw fixed design of pressure gauge assures safety and reliability of the seal.
- The inflator is equipped with fluorescent gauge, which contributes to accuracy of reading.
- Ergonomic handle is comfortable and convenient to grasp.
- The inflator is of good leakproofness since it uses materials with high strength and safety.

Specifications

Model	Nominal Capacity	Maximum Pressure	Syringe Features
BI-20A-10	20cc	30atm	With 10ml Syringe



Add: No.27 Xinke West Road, Luoyang Town, Wujin District,
Changzhou, Jiangsu, China.

Tel: 0086-519-88520118

Fax: 0086-519-88551179-8888

E-mail: sales@czmed.com

Web: <http://www.czmed.com>