



MR. THERMO



**NINGBO HI-LIFE
MEDICAL
TECHNOLOGY CO.,LTD.**

www.hilife-med.com



COMPANY PROFILE

As one of the leading manufacturer in clinical Thermometer industry, Ningbo Hi-life Medical Technology Co.,Ltd, founded in 1992 (Ningbo Huahui Medical Instruments,.Co.,Ltd),and we specialized in developing and manufacturing of medical thermometers,from Mercury thermometers to new Eco-friendly mercury-free thermometers,from digital thermometer to Infrared Forehead thermometer,all of our products are tested strictly to assure safety,accuracy and reliability.

Measuring your body temperature with our professional and Appreciation!





CLINICAL GLASS THERMOMETER

Mercury-Free




6000-A

Large Size, Scale in C, F, C/F



7000-A

Medium Size, Scale in C Only



8000-A

Small Size, Scale in C only



DIGITAL THERMOMETER

Rigid Type



ECT-1

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time (Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



MT-808

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time (Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



MT-811

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time (Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



MT-518

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time (Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



MT-905

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time (Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



MT-909

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time (Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery





DIGITAL THERMOMETER

Cartoon Cap



MT-801C

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time
(Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



Flexible Type



MT-519

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time
(Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



MT-801

Remark:

1. Display range: 32.0~42.9°C
2. Accuracy: $\pm 0.1^{\circ}\text{C}$
3. Min. Scale: 0.1°C
4. Measurement Time
(Reference Only): 1 min $\pm 20\text{sec}$
5. Beeper function
6. Auto Power-off
7. Battery: 1.5V Button Battery



INFRARED THERMOMETER

Forehead / Ear Thermometer



ET-01

Remark:

1. Display range: 32.0~42.9°C
2. Measurement range: 34.0°C -42.9°C
3. Accuracy: 35.0°C -42.0°C / $\pm 0.2^{\circ}\text{C}$
35.0°C below, 42.0°C above / $\pm 0.3^{\circ}\text{C}$
4. Memory function: Memorize the last 9 temperature readings
5. Measurement position: Ears
6. Operating environment: Operating temperature: 5°C -40°C
Relative humidity: $\leq 85\%$
Atmospheric pressure: 700hPa-1060hPa
7. Power: DC 3V(2x1.5V AAA Battery)
8. Service life: 5 years(from manufacturing date)
9. Size: 134X60X37 mm
10. Weight: 60g(battery included)



ET-02

Remark:

1. Display range: 32.0~42.9°C
2. Measurement range: 34.0°C -42.9°C
3. Accuracy: 35.0°C -42.0°C / $\pm 0.2^{\circ}\text{C}$
35.0°C below, 42.0°C above / $\pm 0.3^{\circ}\text{C}$
4. Memory function: Memorize the last 9 temperature readings
5. Measurement position: Ears
6. Operating environment: Operating temperature: 5°C -40°C
Relative humidity: $\leq 85\%$
Atmospheric pressure: 700hPa-1060hPa
7. Power: DC 3V(2x1.5V AAA Battery)
8. Service life: 5 years(from manufacturing date)
9. Size: 134X60X37 mm
10. Weight: 60g(battery included)



INFRARED THERMOMETER

Forehead Thermometer

FHT-1



Model :

- 1.Measuring range: Body temperature: 32.0°C-43.0°C
- 2.Accuracy: 32.0°C-34.9°C: $\pm 0.3^{\circ}\text{C}$
35.0°C-42.0°C: $\pm 0.2^{\circ}\text{C}$
42.1°C-43.0°C: $\pm 0.3^{\circ}\text{C}$
- 3.Display resolution: 0.1°C
- 4.Measuring distance: 1-3cm
- 5.Measurement position: Forehead
- 6.Operating environment: Operating temperature: 15°C-40°C
Relative humidity: $\leq 85\%$
- 7.Atmospheric Pressure: 700hPa-1060hPa
- 8.Storage environment:Storage temperature: - 10°C-55°C
Relative humidity: $\leq 95\%$
Atmospheric Pressure: 700hPa-1060hPa

- 9.Battery: DC 3V (2x1.5V AAA Battery)
- 10.Battery life: About 5000 times can be used
- 11.Memory function: Stores up to 32 groups of measurement values
- 12.Indication: If the ambient temperature is not within 15°C-40°C, the LCDscreen will display "Err"Measuring temperature: >43.0°C, displays "Hi";<32°C displays "Lo"
- 13.LCD screen:Green Light: 32.0°C -37.4°C
Yellow Light: 37.5°C -37.9°C
Red Light: 38.0°C -43.0°C
- 14.Service life:2 years (from manufacturing date)
- 15.Size:138mm×95mm×40mm

RT-101



Model :

- 1.Measuring range: Body temperature: 32.0°C-43.0°C
- 2.Accuracy: 32.0°C-34.9°C: $\pm 0.3^{\circ}\text{C}$
35.0°C-42.0°C: $\pm 0.2^{\circ}\text{C}$
42.1°C-43.0°C: $\pm 0.3^{\circ}\text{C}$
- 3.Display resolution: 0.1°C
- 4.Measuring distance: 1-3cm
- 5.Measurement position: Forehead
- 6.Operating environment: Operating temperature: 15°C-40°C
Relative humidity: $\leq 85\%$
- 7.Atmospheric Pressure: 700hPa-1060hPa
- 8.Storage environment:Storage temperature: - 10°C-55°C
Relative humidity: $\leq 95\%$
Atmospheric Pressure: 700hPa-1060hPa

- 9.Battery: DC 3V (2x1.5V AAA Battery)
- 10.Battery life: About 5000 times can be used
- 11.Memory function: Stores up to 32 groups of measurement values
- 12.Indication: If the ambient temperature is not within 15°C-40°C, the LCD screen will display "Err"Measuring temperature: >43.0°C, displays "Hi";<32°C displays "Lo"
- 13.LCD screen:Green Light: 32.0°C -37.4°C
Yellow Light: 37.5°C -37.9°C
Red Light: 38.0°C -43.0°C
- 14.Service life:2years (from manufacturing date)
- 15.Size:165mm×46mm×47mm

RT-105



Model :

- 1.Measurement range: 32.0° C-43.0° C
- 2.Accuracy: 32.0° C-34.9° C: $\pm 0.3^{\circ}\text{C}$
35.0° C -42.0° C: $\pm 0.2^{\circ}\text{C}$
42.1° C-43.0° C: $\pm 0.3^{\circ}\text{C}$
- 3.Display resolution: 0.1° C
- 4.Measuring distance: 0cm-5cm
- 5.Measurement position: Forehead
- 6.Operating environment: Operating temperature:5° C- 40° C
Relative humidity: $\leq 85\%$
- 7.Atmospheric Pressure: 700hPa- 1060hPa
- 8.Storage environment: Storage temperature:-10° C-55° C
Relative humidity: $\leq 95\%$
- 9.Atmospheric Pressure: 700hPa-1060hPa
- 10.Power: DC3V (2x1.5V AAA Battery)
- 11.Memory function: Stores up to 9 groups of measurement valuesService life:5 years

DIGITAL BATH THERMOMETER



BT-01



BT-02



BT-03



NINGBO HI-LIFE MEDICAL TECHNOLOGY CO.,LTD.

Add: QiaoLi Village, JiangShan Town, Yinzhou District,
Ningbo City, ZheJiang Province, P.R. China 315191

Tel: +86 574 88071300

Fax: +86 574 88451521

E-mail: byron.pan@hilife-med.com, info@hilife-med.com

Http: //www.hilife-med.com